

ERK2 peptide

Cat. No. GTX27992

Applications	IP
Species	Rat

Package
100 µg

Applications

Application Note

Blocking peptide for rabbit polyclonal antibody to ERK2(GTX27948).

Properties

Form	Liquid
Buffer	PBS, BSA
Preservative	0.1% Sodium azide
Storage	Keep as concentrated solution, aliquot and store at 4°C.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Peptide mapping to the carboxy terminus of rat ERK2.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).