

c-Fos peptide

Cat. No. GTX27997

Applications	IP
Species	Human

Package
100 µg

Applications

Application Note

Blocking peptide for rabbit polyclonal antibody to c-Fos(GTX27963).

Properties

Form	Liquid
Buffer	PBS, BSA
Preservative	0.1% Sodium azide
Storage	Keep as concentrated solution, aliquot and store at 4°C.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Peptide mapping to the amino terminus of human c-Fos.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).