

# STAT5a peptide

**Cat. No. GTX28003**

|                     |       |
|---------------------|-------|
| <b>Applications</b> | IP    |
| <b>Species</b>      | Mouse |

**Package**  
100 µg

## Applications

### Application Note

Blocking peptide for rabbit polyclonal antibody to STAT5a(GTX27968).

## Properties

|                        |  |
|------------------------|--|
| <b>Form</b>            | Liquid   |
| <b>Buffer</b>          | PBS, BSA   |
| <b>Preservative</b>    | 0.1% Sodium azide  |
| <b>Storage</b>         | Keep as concentrated solution, aliquot and store at 4°C.                   |
| <b>Concentration</b>   | 0.2 mg/ml (Please refer to the vial label for the specific concentration.) |
| <b>Region/Sequence</b> | Peptide mapping to the carboxy terminus of mouse Stat5a.                   |
| <b>Conjugation</b>     | Unconjugated   |

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).