

p27 Kip1 peptide

Cat. No. GTX28016

Applications	IP
Species	Human

Package
100 µg

Applications

Application Note

Blocking peptide for rabbit polyclonal antibody to p27 Kip1(GTX27961).

Properties

Form	Liquid
Buffer	PBS, BSA
Preservative	0.1% Sodium azide
Storage	Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Peptide mapping to the carboxy terminus of human p27.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).