

Human Albumin protein (Rhodamine)

Cat. No. GTX28031-pro**Species** Human**Package**

1 mg

Properties

Form Lyophilized powder**Buffer** Reconstitute with 1.0ml ddH₂O. Lyophilized from 0.02M Potassium Phosphate, 0.15M NaCl (pH7.2), 10mg/ml BSA (IgG, Protease free).**Preservative** 0.01% Sodium azide**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.**Expression System** Normal Human Serum**Conjugation** Rhodamine [Wavelength](#)
Rhodamine (TRITC)**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).