

Human Albumin protein (HRP)

Cat. No. GTX28032-pro

Species Human

Package
1 mg

Properties

Form Lyophilized powder

Buffer Reconstitute with 1.0ml ddH₂O to 1.0mg/ml. Lyophilized from 0.02M Potassium Phosphate, 0.15M NaCl (pH7.2), 10mg/ml Polyethylene Glycol.

Preservative 0.01% Gentamicin Sulfate

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Expression System Normal Human Serum

Conjugation Horseradish peroxidase(HRP)

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).