

Human Albumin protein (HRP)

Cat. No. GTX28032-pro

Species	Human
----------------	-------

Package	1 mg
----------------	------

PROPERTIES

Form	Lyophilized powder
Buffer	Reconstitute with 1.0ml ddH ₂ O to 1.0mg/ml. Lyophilized from 0.02M Potassium Phosphate, 0.15M NaCl (pH7.2), 10mg/ml Polyethylene Glycol.
Preservative	0.01% Gentamicin Sulfate
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Expression System	Normal Human Serum
Conjugation	Horseradish peroxidase(HRP)
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).