

HLA-DR + HLA-DP antibody [MEM-136]

Cat. No. GTX28094

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
FCM	Assay dependent
IP	Assay dependent

Not tested in other applications.

Product Note This antibody reacts with alpha/beta dimer as well as with dissociated beta-subunit.

Properties


Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	PHA-activated peripheral blood lymphocytes.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX28094 WB Image

WB analysis of Raji cell lysate using GTX28094 HLA-DR + HLA-DP antibody [MEM-136].

Lane 1 : GTX28094

Lane 2 : mouse IgG1 isotype control


GTX28094 IHC-P Image

IHC-P analysis of human small intestine tissue using GTX28094 HLA-DR + HLA-DP antibody [MEM-136].



For full product information, images and publications, please visit our [website](https://www.genetex.com).