

IL1RAP antibody

Cat. No. GTX28109

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5µg/ml
ICC/IF	5µg/ml
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	65 kDa. (Note)
Product Note	No crossreactivity to other members in the IL-1 receptor family.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to 22 amino acids near the carboxy terminus of human Interleukin 1 Receptor AcP.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX28109 WB Image

WB analysis of HeLa cell lysate using GTX28109 IL1RAP antibody.

Dilution : 1µg/ml



GTX28109 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX28109 IL1RAP antibody.

Dilution : 5µg/ml



For full product information, images and publications, please visit our [website](#).