

Survivin peptide

Cat. No. GTX28122

Applications Neutralizing/Inhibition

Species Mouse

Package

50 µg

Applications

Application Note

The peptide is used for blocking the antibody activity of survivin (GTX28113). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C

Properties

Form Liquid

Buffer PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store the peptide at -20°C, stable for one year.

Concentration 0.20 mg/ml (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).