

# Pepsin (ready to use)

**Cat. No. GTX28194**

**Applications** IHC-P

References ( 4 )  
Package  
50 ml

## PRODUCT

**Summary** Antigen retrieval solution-enzymatic (pepsin)

## Applications

### Application Note

This antigen retrieval solution is provided in a ready-to-use format requiring no further preparation. 1.) Apply this solution to tissue sections following their initial hydration step. Incubate for 10 min at 37°C or 45°C (this incubation time may vary depending on the degree of fixation of the tissue sections. 2.) Wash tissue sections well with IHC wash buffer (GTX28195), or PBS or TBS if immunostaining by hand.

## Properties

**Form** Liquid

**Storage** Refridgerate at 4°C. DO NOT FREEZE. Avoid extended exposure to room temperature.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).