

Antigen retrieval condition-Pepsin (ready to use)

Cat No. GTX28194

Application IHC-P, IHC

Reference (4)
 Package
 50 ml

PRODUCT

Summary Antigen retrieval solution-enzymatic (pepsin)

APPLICATION

Application Note

This antigen retrieval solution is provided in a ready-to-use format requiring no further preparation. 1.) Apply this solution to tissue sections following their initial hydration step. Incubate for 10 min at 37°C or 45°C (this incubation time may vary depending on the degree of fixation of the tissue sections. 2.) Wash tissue sections well with IHC wash buffer (GTX28195), or PBS or TBS if immunostaining by hand.

PROPERTIES

Form Liquid

Storage Refrigerate at 4°C. DO NOT FREEZE. Avoid extended exposure to room temperature.

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).