

Adenovirus hexon antibody [8C4]

Cat. No. GTX28249

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, Immunodiffusion, IHC, Lateral Flow, Sandwich ELISA
Reactivity	Adenovirus

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Immunodiffusion	Assay dependent
IHC	Assay dependent
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent
Note : Capture : GTX28249, Detection : GTX28251	

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Hexon antigen of human adenovirus
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).