

C Reactive Protein antibody [C5]

Cat. No. GTX28279

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Application	ELISA, IHC, Purification, Turbidimetry, Sandwich ELISA
Reactivity	Human

Package
100 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
IHC	Assay dependent
Purification	Assay dependent
Turbidimetry	Assay dependent
Sandwich ELISA	Assay dependent
Note : Capture : GTX28279, Detection : GTX28278	

Not tested in other applications.

Calculated MW 25 kDa. ([Note](#))

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human CRP, derived from pleural/ascetic fluid or plasma.
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 02 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).