

# FK506 antibody [FK1]

**Cat. No. GTX28316**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgM
<b>Applications</b>	ELISA
<b>Reactivity</b>	Species independent

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

**Product Note** This antibody recognizes conjugated and non-conjugated FK-506 with apparent dissociation constant  $K_d = 10^{-9}$  M.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	50mM Tris pH8.0, 150mM NaCl
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Conjugate FK-506-BSA.
<b>Purification</b>	Purified by euglobulin precipitation and dialysis
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).