

FK506 antibody [FK1]

Cat. No. GTX28316

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	ELISA
Reactivity	Species independent

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note This antibody recognizes conjugated and non-conjugated FK-506 with apparent dissociation constant $K_d = 10^{-9}$ M.

Properties

Form	Liquid
Buffer	50mM Tris pH8.0, 150mM NaCl
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Conjugate FK-506-BSA.
Purification	Purified by euglobulin precipitation and dialysis
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).