

Interferon gamma antibody [H3-1]

Cat. No. GTX28318

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

This antibody can be used to neutralize the anti-viral activity of interferon gamma. Can also be used for the isolation and purification of recombinant interferon gamma.

Product Note Reacts with human interferon gamma.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE. We recommend to parafilm vials during storage to prevent evaporation of buffer.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length protein.
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).