

RAIDD antibody

Cat. No. GTX28426

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1µg/ml
ICC/IF	5µg/ml
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide corresponding to aa 99- 117 of human RAIDD (accession no. AAB42217).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx28426 ICC/IF Image

Immunocytochemical staining of RAIDD in HeLa cells with RAIDD antibody (GTx28426) at 5µg/ml.



GTx28426 WB Image

Western blot analysis of RAIDD in whole cell lysates from HeLa (H) or K562 (K) cells with Anti-RAIDD (IN) (Cat. No. GTx28426) at 1:500 dilution.



GTx28426 ICC/IF Image

ICC/IF analysis of HeLa cells using GTx28426 RAIDD antibody.
Dilution : 20 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).