

MMP2 antibody

Cat. No. GTX28481

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent

Not tested in other applications.

Calculated MW 74 kDa. ([Note](#))

Product Note Identity has been confirmed by Western blot and by double diffusion (Ouchterlony) against human MMP2 and a known anti-humanMMP2.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	ProMMP-2 purified from human cultured rheumatoid synovial fibroblasts, shown to be homogeneous by SDS-PAGE.
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).