

Human Lung: Left Lower Lobe (Normal) tissue lysate, Cytoplasmic Fraction

Cat. No. GTX28691

Species Human

Package
100 µg

PRODUCT

Summary

Human lung left lower lobe tissue cytoplasmic protein lysate was prepared by isolating the cytoplasmic protein from whole tissue homogenates using a proprietary technique. The human lung left lower lobe tissue was frozen in liquid nitrogen immediately after excision and then stored at -70°C. For quality control purposes, the isolated lung left lower lobe tissue cytoplasmic protein pattern on SDS-PAGE gel is shown to be consistent for each lot by visualization with coomassie blue staining.

Tissue type Normal

Properties

Form Liquid

Buffer HEPES, MgCl₂, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).