

# Human Brain: Thalamus (Normal) tissue lysate

**Cat. No. GTX28710**

**Species** Human

**References** ( 1 )

**Package**

200 µg

## PRODUCT

**Summary** Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.

**Tissue type** Normal

## Properties

**Form** Liquid

**Buffer** HEPES, MgCl<sub>2</sub>, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration** 5 mg/ml (Please refer to the vial label for the specific concentration.)

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).