

Ceruloplasmin antibody

Cat. No. GTX28749

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	IP
Reactivity	Human

Package
500 µl

Applications

Application Note

Applications for immunohistochemistry have not been determined

Product Note This antibody reacts with human ceruloplasmin, a 140 kD a globulin of blood plasma.

Properties

Form	Liquid
Buffer	20mM Sodium Phosphate, 140mM NaCl
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Full length protein (Human).
Purification	Whole antiserum
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).