

TFF1 antibody

Cat. No. GTX28762

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr, ELISA
Reactivity	Human

References (1)

Package

250 µl

Applications

Application Note

Immunohistochemistry 1:25-1:50. This antibody can be used on both frozen and formalin-fixed, paraffin-embedded tissue sections. Prolonged fixation in buffered formalin can destroy the epitope. This antibody may be used at a dilution of 1:25-1:50. For paraffin sections, incubation of primary antibody for one hour at 37° is recommended. The optimum conditions should be determined by the individual laboratory.

Product Note This antibody is specific to pS2 protein, also known as pNR-2.

Properties

Form	Liquid
Buffer	Serum
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Synthetic peptide representing a highly antigenic site near the C-terminus of the pS2 protein.
Purification	Whole antiserum
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).