

Jurkat cell lysate, Membrane Fraction

Cat. No. GTX28894
Species Human

Package
100 µg

PRODUCT

Summary

Jurkat (Human Acute T cell Leukemia) cell membrane protein lysate was prepared by isolating the membrane protein from whole tissue homogenates using a proprietary technique. The Jurkat (Human Acute T cell Leukemia) cell was frozen in liquid nitrogen immediately after excision and then stored at -70°C. For quality control purposes, the isolated Jurkat (Human Acute T cell Leukemia) cell membrane protein pattern on SDS-PAGE gel is shown to be consistent for each lot by visualization with coomassie blue staining.

Properties

Form Liquid

Buffer HEPES, MgCl₂, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).