

Placental Alkaline Phosphatase antibody

Cat. No. GTX28942

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA, EM
Reactivity	Human

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
EM	Assay dependent

Not tested in other applications.

Product Note

This antibody recognizes alkaline phosphatase, it has been found to be reactive with both placental and bone alkaline phosphatase.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified alkaline phosphatase from human placenta.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).