

Human WBC (Normal) lysate

Cat. No. GTX29043

Species Human

Package
100 µg

PRODUCT

Summary Peripheral Blood Leukocyte (PBL's) tissue lysate was prepared by homogenization in lysis buffer. Tissue and cell debris was removed by centrifugation. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Tissue type Normal

Properties

Form Liquid

Buffer 1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 2 mg/ml (Please refer to the vial label for the specific concentration.)

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).