

Angiopoietin 2 peptide

Cat. No. GTX29077

Applications	Neutralizing/Inhibition
Species	Mouse

Package
50 µg

Applications

Application Note

For BL: Control peptides should be used at 1.0 µg per 1.0 µl of anti-serum in per assay. Optimal dilutions/concentrations should be determined by the researcher.

Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purity	> 95%
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).