

Human Pancreas (Tumor) tissue lysate

Cat. No. GTX29129

Species	Human
---------	-------

Package	100 µg
---------	--------

PRODUCT

Summary

Pancreas tissue lysate was prepared by homogenization in modified RIPA buffer. Tissue and cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The product was boiled for 5 min in 1 x SDS sample buffer containing 5% β-mercaptoethanol.

Tissue type

Tumor

Properties

Form

Liquid

Buffer

1x SDS sample buffer (50mM Tris-HCl, 12.5% Glycerol, 1% Sodium Dodecylsulfate, 0.01% Bromophenol blue), 5% β-mercaptoethanol

Preservative

No preservatives

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration

2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 1