

Goat Anti-Human IgE antibody

Cat. No. GTX29159

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, ELISA, IHC
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	10 - 2 ug/ml
ELISA	Assay dependent
IHC	10 -2 ug/ml

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human IgE myeloma proteins.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX29159 ELISA Image

ELISA analysis of immunogen using GTX29159 Human IgE antibody.


GTX29159 Image

SDS-PAGE of GTX29159 Human IgE antibody.

Lane 1 : Primary antibody

Lane 2 : Marker

Loading : 0.2μg



For full product information, images and publications, please visit our [website](https://www.genetex.com).