

Goat Anti-Mouse IgG2c antibody

Cat. No. GTX29168

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, ELISA
Reactivity	Mouse

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 0.05 ug/ml
IP	2 ug/mg
ELISA	2 ug/ml

Not tested in other applications.

Product Note

The affinity purified antibodies react with IgG2c in C57BL/6, SJL, C57BL/10, CB20, C57BL/6 by Balb/C crosses and pools of serum of outbred mice obtained from several commercial sources. The anti IgG2c does not react with L-chains or IgG2a, IgG2b, IgG1 or IgG3. This reagent may also be used where a reagent specific for Igh 1-b allele is indicated.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Goats were immunized with highly purified IgG2c isolated from C57BL/6 mice, in FreundOs adjuvant.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX29168 ELISA Image**

ELISA analysis of immunogen using GTX29168 Mouse IgG2c antibody.



For full product information, images and publications, please visit our [website](#).