

Plasmodium falciparum antibody

Cat. No. GTX29206

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	WB, ELISA
Reactivity	Plasmodium falciparum

Package
500 µg

Applications
Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Histidine rich protein II from Plasmodium falciparum
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX29206 ELISA Image

ELISA analysis of immunogen using GTX29206 Plasmodium falciparum antibody.



For full product information, images and publications, please visit our [website](#).