

Triiodothyronine (T3) antibody [ME-124]

Cat. No. GTX29237

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ELISA
Reactivity	Species independent

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note It shows less than 1% cross-reactivity with T2, rT3, and T4. Affinity Constant: 5×10^8

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.60 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	3,5,3'-Triiodothyronine coupled to BSA
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).