

## Human Colon descending (Normal) tissue lysate, Membrane Fraction

Cat. No. GTX29412

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

|         |        |
|---------|--------|
| Package | 100 µg |
|---------|--------|

## PRODUCT

## Summary

Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.

## Tissue type

Normal

## Properties

## Form

Liquid

## Buffer

HEPES, MgCl<sub>2</sub>, KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors

## Preservative

No preservatives

## Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

## Concentration

2 mg/ml (Please refer to the vial label for the specific concentration.)

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).