

PAG antibody [MEM-255]

Cat. No. GTX29527

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-P, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Product Note

This antibody recognizes an epitope (aa 235-280) of Csk-binding protein (Cbp) located in the cytoplasmic domain. We do not recommend use of this product for Mouse,Rat,Bovine samples.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant intracellular fragment (aa 97-432) of human Cbp (PAG).
Purification	Protein A purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx29527 IHC-P Image

IHC-P analysis of germinal centers of lymph node lymphoid follicle and in follicular lymphoma (it is absent from mantle zone) using GTx29527 PAG antibody [MEM-255]. Cbp/PAG is also expressed more weakly in T cells in tonsil and the thymic medulla.



For full product information, images and publications, please visit our [website](#).