

Human MCP4 protein

Cat. No. GTX29560-pro

Species Human

Package
10 µg

PRODUCT

Summary

Monocyte Chemotactic Protein-4 Human Recombinant produced in E.Coli is a non-glycosylated, Polypeptide chain containing 75 amino acids and having a molecular mass of 8.6 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 20mM PBS (pH7.4), 130mM NaCl.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence QPDALNVPSTCCFTFSSKKISLQRLKSYVITTSRCPQKAVIFRTKLGKEICADPKEKWWQNYMKHLGRKAHTLKT.

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).