

# Human IL15 protein

**Cat. No. GTX29583-pro**
**Species** Human

**Package**

10 µg

## PRODUCT

**Summary**

Interleukin-15 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 114 amino acids and having a molecular mass of 12.9 kDa.

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with ddH<sub>2</sub>O but not less than 100µg/ml. Lyophilized from PBS (pH7.4).

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** MNWVNVISDL KKIEDLIQSM HIDATLYTES DVHPSCKVTA MKCFLELQV ISLESGDASI HDTVENLIIL ANNSLSSNGN VTESGCKECE ELEEKNIKEF LQSFVHIVQM FINTS.

**Expression System** E. coli

**Purity** > 98% by SDS-PAGE and HPLC.

**Endotoxin** < 1 EU/µg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).