

Human KGF protein

Cat. No. GTX29601-pro

Species Human

Package

10 µg

PRODUCT

Summary

KeRatinocyte Growth Factor-1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 164 amino acids and having a molecular mass of 18.995 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 20mM PBS (pH8.0), 1M NaCl.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MCNDMTPEQM ATNVNCSSPE RHTRSYDME GGDIVRRRLF CRTQWYLRI KRGKVKGTQE MKNNYNIMEI RTVAVGIVAI KGVSEFYLA MNKEGKLYAK KECNEDCNFK ELLENHYNT YASAKWTHNG GEMFVALNQK GIPVRGKGTK KEQKTAHFLP MAIT.

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).