

Human IL3 protein

Cat. No. GTX29621-pro

Species Human

Package

10 µg

PRODUCT

Summary

Interleukin-3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 133 amino acids and having a molecular mass of 15 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100µg/ml. Lyophilized from 0.3mg/ml NaHCO₃ solution.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Pro-Met-Thr-Gln.

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).