

# Human MCP2 protein

**Cat. No. GTX29672-pro**

**Species** Human

**Package**  
10 µg

## PRODUCT

### Summary

Monocyte Chemotactic Protein-2 Human Recombinant produced in E.Coli is a non-glycosylated, Polypeptide chain containing 76 amino acids and having a molecular mass of 8.904 kDa.

## Properties

**Form** Lyophilized powder

**Buffer** Reconstitute with ddH<sub>2</sub>O but not less than 100µg/ml. Lyophilized from 1.0mg/ml solution with no additives.

**Preservative** No preservatives

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** The sequence of the first five N-terminal amino acids was determined and was found to be Gln-Pro-Asp-Ser-Val.

**Expression System** E. coli

**Purity** > 98% by SDS-PAGE and HPLC.

**Endotoxin** < 1 EU/µg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).