

Human MIP1 Beta protein

Cat. No. GTX29676-pro

Species Human

Package

10 µg

PRODUCT

Summary

Macrophage Inflammatory Protein-1 beta Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 69 amino acids and having a molecular mass of 7.62 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100µg/ml. Lyophilized from 1.0mg/ml solution with no additives.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The sequence of the first five N-terminal amino acids was determined and was found to be, Ala-Pro-Met-Gly-Ser.

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).