

Human RANTES protein

Cat. No. GTX29680-pro

Species Human

Package

20 µg

PRODUCT

Summary

Rantes Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 68 amino acids and having a molecular mass of 7.809 kDa.

Properties

Form	Lyophilized powder
Buffer	Reconstitute with ddH ₂ O but not less than 100 µg/ml. Lyophilized from 0.2µm-filtered 50mM NaAc (pH5.0), 100mM NaCl.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	SPYSSDTTPC CFAYIARLP RAHIKEYFYT SGKCSNPAVV FVTRKNRQVC ANPEKKWVRE YINSLEMS.
Expression System	E. coli
Purity	> 98% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).