

# Human PDGF-AA protein

**Cat. No. GTX29702-pro**

**Species** Human

**Package**  
10 µg

## PRODUCT

**Summary** Platelet-Derived Growth Factor AA Human Recombinant is a homodimeric, non-glycosylated, polypeptide chain containing 2 x 125 amino acids and having a total molecular mass of 28.511 kDa.

## Properties

**Form** Lyophilized powder

**Buffer** Batch dependent (please contact us for details)

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

**Region/Sequence** SIEEAVPAVC KTRTVIYEIP RSQVDPTSAN FLIWPPCDEV KRCTGCCNTS SVKCQPSRVH HRSVKVAKVE YVRKKPKLKE VQVRLEEHL E CACATSLNP DYREEDTGRP RESGKKRK RK RLKPT.

**Expression System** E. coli

**Purity** > 98% by SDS-PAGE and HPLC.

**Endotoxin** < 1 EU/µg

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](https://www.genetex.com).