

Human Brain: Amygdala (Alzheimer's Disease) tissue lysate, Membrane Fraction

Cat. No. GTX29712

Species	Human
---------	-------

Package	100 µg
---------	--------

PRODUCT

Summary	Tissue total protein is prepared from whole tissue homogenates. It presents a consistent pattern on SDS-PAGE analysis. The total protein is stored in a buffer with protease inhibitor cocktail.
---------	--

Tissue type	Neurological Disease
-------------	----------------------

Properties

Form	Liquid
------	--------

Buffer	HEPES, MgCl ₂ , KCl, EDTA, Sodium deoxycholate, Sucrose, NP-40, Glycerol, a cocktail of protease inhibitors
--------	--

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---

Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 16 Page 1 of 1