

Human MIG protein

Cat. No. GTX29721-pro

Species Human

Package

20 µg

PRODUCT

Summary

MIG (monokine induced by gamma-interferon) Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 103 amino acids and having a molecular mass of 11.7 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 20mM PBS (pH7.4), 50mM NaCl.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence TPVVRKGRCSISTNQGTIHLQSLKDLKQFAPSPSCEKIEIATLKNGVQ
TCLNPDSADVKELIKKEKQVSQKKKQKNGKKHKKKVLKVRKSQRSRQKKTT.

Expression System E. coli

Purity > 98% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).