

Human CCL21 protein

Cat. No. GTX29852-pro

Applications	Functional Assay
Species	Human

Package
20 µg

PRODUCT

Summary Exodus-2 Human Recombinant produced in E.Coli as a single, non-glycosylated, polypeptide chain containing 111 amino acids and having a molecular mass of 12.219 kDa.

Properties

Form	Lyophilized powder
Buffer	Reconstitute with ddH ₂ O to 0.1-1.0mg/ml. Lyophilized from 20mM Phosphate buffer (pH7.4), 150mM NaCl.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	SDGGAQDCCL KYSQRKIPAK VRSYRKQEP SLGCSIPAIL FLPRKRSQAE LCADPKELWV QQLMQHLDKT PSPQKPAQGC RKDRGASKTG KKGKGSKGCR KTERSQTPKG P.
Expression System	E. coli
Purity	> 98% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).