

Mouse Interferon Gamma protein

Cat. No. GTX29922-pro

Species Mouse

Package
100 µg

PRODUCT

Summary

Interferon-gamma Mouse Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 134 amino acids and having a molecular mass of 15.6 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O or 20mM AcOH to 0.1-0.5mg/ml. Lyophilized from 0.2µm-filtered PBS (pH7.4), 5% Trehalose.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MHGTVIESLE SLNNYFNSSG IDVEEKSLFL DIWRNWQKDG DMKILQSQII SFYLRLFEVL KDNQAISNNI SVIESHLITT FFSNSKAKKD AFMSIAKFEV NNPQVQRQAF NELIRVVHQL LPESLRKRK RSRC.

Expression System E. coli

Purity > 95% by SDS-PAGE and HPLC.

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).