

Human IL12B / IL12 p40 protein, His tag

Cat. No. GTX29993-pro

Species Human

Package
10 µg

PRODUCT

Summary Interleukin-12 p40 His Human Recombinant produced in E.Coli is single, a non-glycosylated, polypeptide chain containing 306 amino acids fragment (23-328) with an amino-terminal hexahistidine tag.

Properties

Form Liquid

Buffer PBS, 50% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 0.0746 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminus His-Tag

Expression System E. coli

Purity > 95% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 11 Page 1 of 1