

FGF10 antibody

Cat. No. GTX30007

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

References (1)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.2µg/ml
IHC-P	Assay dependent
ELISA	0.5-2µg/ml

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human FGF 10
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX30007 IHC-P Image

IHC-P analysis of human carcinoid tissue using GTX30007 FGF10 antibody. Antigen retrieval using citrate buffer at pH6.0.


GTX30007 IHC-P Image

IHC-P analysis of human carcinoid tissue using GTX30007 FGF10 antibody. Antigen retrieval using citrate buffer at pH6.0.


GTX30007 IHC-P Image

IHC-P analysis of human carcinoid tissue using GTX30007 FGF10 antibody. Antigen retrieval using citrate buffer at pH6.0.



For full product information, images and publications, please visit our [website](https://www.genetex.com).