

Golgin 97 antibody [CDF4]

Cat. No. GTX30100

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, IP, IHC
Reactivity	Human, Mouse, Rat, Dog, Primate

References (2)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1.0 µg/ml
IHC-P	Assay dependent
IP	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 88 kDa. ([Note](#))

Properties

Form	Liquid
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human golgin-97.
Purification	Purified IgG
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx30100 ICC/IF Image

ICC/IF analysis of permeabilized HeLa cells using GTx30100 Golgin 97 antibody [CDF4].

Green : Primary antibody

Blue : cell nucleus



GTx30100 ICC/IF Image

ICC/IF analysis of muntjac skin cells using GTx30100 Golgin 97 antibody [CDF4].

Green : Primary antibody

Purple : Actin

Blue : cell nucleus



GTx30100 ICC/IF Image

ICC/IF analysis of permeabilized HDFn cells using GTx30100 Golgin 97 antibody [CDF4].

Green : Primary antibody

Red : Actin

Blue : cell nucleus



GTx30100 ICC/IF Image

ICC/IF analysis of permeabilized HDFn cells using GTx30100 Golgin 97 antibody [CDF4].

Green : Primary antibody

Red : Actin

Blue : cell nucleus


GTx30100 ICC/IF Image

ICC/IF analysis of permeabilized HeLa cells using GTx30100 Golgin 97 antibody [CDF4].

Green : Primary antibody

Red : Actin

Blue : cell nucleus

Dilution : 1 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).