

HSDL2 antibody [2A11BG6AH6]

Cat. No. GTX30325

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, IP, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	10 µg/ml
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	HEPES buffered saline
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human HSDL2
Purification	Purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx30325 ICC/IF Image

ICC/IF analysis of HDFn cells using GTx30325 HSDL2 antibody [2A11BG6AH6].

Fixation : 4% paraformaldehyde 20 mins

Permeabilization : Triton X-100 (0.1%, 15 min)

Dilution : 10 µg/ml overnight at 4°C



For full product information, images and publications, please visit our [website](#).